



Nutritional Delivery Product Catalogue

Kangaroo™: The complete enteral feeding solution

Quality products from the leader in innovative enteral feeding pump technology



Nutritional Delivery Products

The Kangaroo™ portfolio offers a safe, innovative and complete line of enteral feeding products for use across the continuum of care. The Kangaroo™ brand has a long history of providing reliable and trusted products for clinicians and caregivers to optimize clinical efficiency. Our high-performing products are preferred by both clinicians and procurement professionals.†

Transitioning to ENFit™? Be sure to note the product description throughout the catalog for ENFit™ design. Where applicable, the ENFit™ version of the code is listed next to the non-ENFit™ version of the code for ease of reference during a conversion.

American Society of Parenteral and Enteral Nutrition (ASPEN) 2016 Clinical Guidelines

ASPEN developed standards of practice for nutrition support. The clinical guidelines include:

- Nutrition support therapy in the form of early enteral nutrition (EN) should be initiated within 24 to 48 hours in the critically ill patient who is unable to maintain volitional intake.
- Enteral feeding protocols should be designed and implemented to increase the overall percentage of goal calories provided.
- Feeding tubes should be flushed with water after each feeding and medication administration, in between medications and every four hours during a continuous infusion.

†Customer Driven Growth Survey; March 2016:96-98.

Enteral Feeding Pumps and Pump Sets	4
Kangaroo™ ePump, Joey and Connect Enteral Feeding Pumps	4
Extension Sets, Gravity Sets and Irrigation Kits	9
Kangaroo™ Extension Sets	9
Kangaroo™ Gravity Sets	9
Kangaroo™ I.V. Pouch and Irrigation Kits	9
Enteral Access Devices	10
Kangaroo™ IRIS Technology	10
Kangaroo™ Nasogastric Feeding Tubes	11
Kangaroo™ Naso-Jejunal Feeding Tubes	12
Kangaroo™ Neonatal and Pediatric Feeding Tubes	13
Kangaroo™ Neonatal and Pediatric Feeding System	14
Kangaroo™ Gastrostomy Tubes	16
Kangaroo™ Safety P.E.G. Kits	17
Kangaroo™ Skin Level Tubes	18
Accessories	19
Pharmacy Solutions	20
Monoject™ Oral/Enteral Syringes	20
Monoject™ Syringe Tip Caps	21
Monoject™ Medicine Straws	21
Monoject™ UV/Light Inhibiting Transport Bags	21
Monoject™ ENFit™ Universal Bottle Adapters	21
GI Products	22
Salem Sump™ Tubes	22
Stomach Tubes	24
CO ₂ Detectors	26
Index	28

Kangaroo™ ePump, Joey and Connect

Kangaroo™ ePump Enteral Feeding Pump

Features feed and flush technology that automates flushing for tube patency and delivery of hydration.

- Fully programmable variable feed and flush pump
- Flow rate: 1 to 400 mL/hr in 1 mL increments
- Flush rate: 1960 mL/hr
- Flush volume: 10 to 500 mL/hr in 1 mL increments
- Flush interval: 1 to 24 hr in 1 hr increments
- Multilingual: 15 languages
- Battery life: 15 hours
- Weight without pole clamp: 2.4 lb.
- Weight with pole clamp: 2.7 lb.
- Feeding history up to 72 hours
- Attitude independent
- Intuitive user interface
- Auto-prime



Z482400J



Z773656



Z773662



Z716154

Item No.	Description	Packaging
Pump		
Z482400J	Kangaroo™ ePump with Pole Clamp	1 ea/ca
Pump Accessories		
Z382491	Kangaroo™ ePump Power Adapter	1 ea/ca
Z382492	Kangaroo™ ePump Pole Clamp	1 ea/ca
Z1054715	Kangaroo™ ePump Battery	1 ea/ca

Kangaroo™ ePump Enteral Feeding Sets

- Anti-free flow design
- DEHP-free
- ENFit™ connection

Item No.	Description	Packaging
Pump sets		
Z772055	Kangaroo™ ePump 500 mL Pump Set	30 ea/ca
Z773656	Kangaroo™ ePump 1000 mL Pump Set	30 ea/ca
Z674668	Kangaroo™ ePump 500 mL Feed Set with 500 mL Flush Bag	30 ea/ca
Z773662	Kangaroo™ ePump 1000 mL Feed Set with 1000 mL Flush Bag	30 ea/ca
Z716154	Kangaroo™ ePump Sterile 100 mL Burette Set	30 ea/ca
Z775659	Kangaroo™ ePump ENPlus Spike Set	30 ea/ca
Z775100	Kangaroo™ ePump ENPlus Spike with 1000 mL Flush Bag	30 ea/ca
Z776150	Kangaroo™ ePump 100 mL Burette Recertification Set (not ENFit™)	5 ea/ca

Kangaroo™ Joey Enteral Feeding Pump

Features feed and flush technology that automates flushing for tube patency and delivery of hydration with a compact, portable design.

- Fully programmable variable feed and flush pump
- Feed rate: 1 to 400 mL/hr in 1 mL increments
- Flush volume: 10 to 500 mL/hr in 1 mL increments
- Flush interval: 1 to 24 hr in 1 hr increments
- Feeding history up to 72 hours
- Battery life: 18 hours
- Weight without pole clamp: 1.7 lb.
- Weight with pole clamp: 2.3 lb.
- Multilingual: 15 languages
- Ambulatory-specific feed set design
- Tubing length 9 in. (from the bottom of the bag to the distal end of the set)
- Attitude independent
- Intuitive user interface
- Auto-prime



Z383400

Item No.	Description	Packaging
Pump		
Z383400	Kangaroo™ Joey Enteral Feeding Pump with Pole Clamp	1 ea/ca
Pump Accessories		
Z383493	Kangaroo™ Joey Pole Clamp	1 ea/ca
Z383491	Kangaroo™ Joey Power Cord	1 ea/ca
ZF010506	Kangaroo™ Joey Battery	1 ea/ca



Z763656

Kangaroo™ Joey Enteral Feeding Pump Sets

- Anti-free flow design
- DEHP-free
- ENFit™ connection

Item No.	Description	Packaging
Pump sets		
Z762055	Kangaroo™ Joey 500 mL Pump Set	30 ea/ca
Z763656	Kangaroo™ Joey 1000 mL Pump Set	30 ea/ca
Z763662	Kangaroo™ Joey 1000 mL Pump Set with Flush Bag	30 ea/ca
Z765100	Kangaroo™ Joey Safety Screw Spike with 1000 mL Flush Bag	30 ea/ca
Z765559	Kangaroo™ Joey ENPlus Spike Set	30 ea/ca



Z763662



Z765559

Kangaroo™ Joey Backpacks

Kangaroo™ Joey Backpacks are specially designed and tested to hold Kangaroo™ Joey Enteral Feeding Pumps.

Item No.	Description	Packaging
Backpacks		
Z770026	Kangaroo™ Joey Blue Mini Backpack	1 ea/ca
Z770028	Kangaroo™ Joey Blue Super Mini Backpack	1 ea/ca
Z770034	Kangaroo™ Joey Pink Mini Backpack	1 ea/ca
Z770029	Kangaroo™ Joey Pink Super Mini Backpack	1 ea/ca
Z770027	Kangaroo™ Joey Green Mini Backpack	1 ea/ca
Z770032	Kangaroo™ Joey Green Super Mini Backpack	1 ea/ca
Z770025	Kangaroo™ Joey Black Mini Backpack	1 ea/ca
Z770031	Kangaroo™ Joey Black Super Mini Backpack	1 ea/ca



Z384400

Kangaroo™ Connect Enteral Feeding Pump and Accessories

Designed for the ambulatory patient, the Kangaroo™ Connect Enteral Feeding Pump is the first wirelessly connected enteral feeding pump on the market and is the smallest, lightest and quietest Kangaroo™ pump available.

- Most accurate Kangaroo™ enteral feeding pump at +/- 5% accuracy
- Expanded feed rate range
- Intuitive color screen with clear visual notifications
- Up to 24-hour battery life
- 72 hours of feeding history
- Travel accessories including car charger
- Water resistance rating of IPX6 enables easier cleaning
- Simplified recertification
- Auto-prime

Item No.	Description	Packaging
Z384400	Kangaroo™ Connect Enteral Feeding Pump with Pole Clamp and Power Adapter	1 ea/ca
Z384491	Kangaroo™ Connect Power Cord with Adapter	1 ea/ca
Z384492	Kangaroo™ Connect Pole Clamp	1 ea/ca
Z384494	Kangaroo™ Connect Car Charger	1 ea/ca

Kangaroo™ Connect Enteral Feeding Sets for Regular Pump — Not for Use with Wireless Connectivity

- Anti-free flow design
- DEHP-free
- ENFit™ connection
- Cassette pump set design with auto prime feature

Item No.	Description	Packaging
Z77000FD	Kangaroo™ Connect ENPlus Spike Set, Non-sterile	30 ea/ca
Z77100FD	Kangaroo™ Connect 1000 mL Bag Set, Non-sterile	30 ea/ca
Z77500FD	Kangaroo™ Connect 500 mL Bag Set, Non-sterile	30 ea/ca



Z77100FD

Kangaroo™ Connect Enteral Feeding Sets — For Use with Wireless Connectivity

Item No.	Description	Packaging
Z1KCP	Kangaroo™ Connect ENPlus Spike ENTelliSet with ENFit™ Connector, Non-sterile	30 ea/ca
Z1000KCP	Kangaroo™ Connect ENTelliset with ENFit™ Connector, 1000 mL, Non-sterile	30 ea/ca
Z500KCP	Kangaroo™ Connect ENTelliset with ENFit™ Connector, 500 mL, Non-sterile	30 ea/ca



Z1000KCP

Kangaroo™ Connect Hubs for Wireless Connectivity

Item No.	Description	Packaging
Z384500CBW	Kangaroo™ Connect WiFi/Ethernet Communications Hub	1 ea/ca
Z384500CBC	Kangaroo™ Connect Cellular Communications Hub	1 ea/ca



Z384500CBW



Z384500CBC

Kangaroo™ Connect Backpacks



Z770035M shown with Z384400

Item No.	Description	Packaging
Z770035L	Kangaroo™ Connect Backpack, Purple, Large	1 ea/ca
Z770037L	Kangaroo™ Connect Backpack, Black, Large	1 ea/ca
Z770035M	Kangaroo™ Connect Backpack, Purple, Medium	1 ea/ca
Z770037M	Kangaroo™ Connect Backpack, Black, Medium	1 ea/ca
Z770035S	Kangaroo™ Connect Backpack, Purple, Small	1 ea/ca
Z770037S	Kangaroo™ Connect Backpack, Black, Small	1 ea/ca



Extension Sets, Gravity Sets and Irrigation Kits

Kangaroo™ Extension Sets

- Non-sterile

Non-ENFit™ code	ENFit™ code	Description	Packaging
Z8884705008	Z8884705008E	Y-site Extension Set	30 ea/ca
Z8884700207	Z8884700207E	Pump Set 4 ft. Tubing Extension	30 ea/ca

Kangaroo™ Gravity Sets

- Non-sterile
- ENFit™ connection

Item No.	Description	Packaging
Z702505	1000 mL Large Bore Gravity Set	30 ea/ca
Z8884702500	1000 mL Gravity Bag	30 ea/ca



Z702505

Kangaroo™ ENFit™ Transition Connector

- Non-sterile

Item No.	Description	Packaging
Z765000	Kangaroo™ ENFit™ Transition Connector	20 ea/ca

Kangaroo™ I.V. Pouch and Irrigation Kits

- Convenient tear-to-open packaging, resealable flap
- Transparent, polyethylene bag hangs easily, clearly labeled
- Non-sterile

Non-ENFit™ code	ENFit™ code	Description	Packaging
Z84065	Z84065E	I.V. Pouch with 60 cc Piston Syringe	30 ea/ca
Z8884700108	Z8884700108E	Irrigation Kit Syringe with 500 cc Container and Basin	20 ea/ca
Z8884700116	Z8884700116E	Irrigation Kit Syringe with 500 cc Container	30 ea/ca




Z8884700108E

IRIS Technology

Kangaroo™ Feeding Tube with IRIS Technology

The first disposable feeding tube with an integrated real-time imaging system (IRIS) to visually aid in the placement of small-bore feeding tubes.

- Features a camera integrated with a small-bore feeding tube to allow clinicians to identify anatomical markers during the placement of the tube



Z386100

Item No.	Description	Packaging
Console		
Z386100	Kangaroo™ IRIS Console with Accessories (Global Power Adapter, Battery, Interface Cable and Pole Clamp)	1 ea/ca
Feeding Tubes		
Z460836E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 36 in. (91 cm)	5 ea/ca
Z460843E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 43 in. (109 cm)	5 ea/ca
Z460855E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 55 in. (140 cm)	5 ea/ca
Z461036E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 36 in. (91 cm)	5 ea/ca
Z461043E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 43 in. (109 cm)	5 ea/ca
Z461055E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 55 in. (140 cm)	5 ea/ca
Z461236E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 36 in. (91 cm)	5 ea/ca
Z461243E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 43 in. (109 cm)	5 ea/ca
Z461255E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 55 in. (140 cm)	5 ea/ca
Accessories		
Z388101	Kangaroo™ IRIS Power Adapter	1 ea/ca
Z388102	Kangaroo™ IRIS Replacement Battery	1 ea/ca
Z388103	Kangaroo™ IRIS Interface Cable	1 ea/ca
Z388104	Kangaroo™ IRIS Mounting Clamp	1 ea/ca
Z388105E	Kangaroo™ IRIS Insufflation Device	5 ea/ca

Nasogastric Feeding Tubes

Kangaroo™ Feeding Tubes with Dobbhoff™ Tip

- Weighted feeding tubes with and without flow-through stylets
- Radiopaque polyurethane material
- Centimeter markings
- MRI conditional once stylet is removed
- DEHP-free
- ENFit™ connection

Non-ENFit™ code	ENFit™ code	Includes stylet	Weight	French size	Length	Packaging
Z8884710842	Z8884710842E	No	7 g	8 Fr (2.7 mm)	43 in. (109 cm)	10 ea/ca
Z8884710859	Z8884710859E	Yes	7 g	8 Fr (2.7 mm)	43 in. (109 cm)	10 ea/ca
Z8884710826	Z8884710826E	Yes	7 g	8 Fr (2.7 mm)	55 in. (140 cm)	10 ea/ca
Z8884711006	Z8884711006E	Yes	7 g	10 Fr (3.3 mm)	43 in. (109 cm)	10 ea/ca
Z8884711246	Z8884711246E	No	7 g	12 Fr (4.0 mm)	43 in. (109 cm)	10 ea/ca
Z8884711253	Z8884711253E	Yes	7 g	12 Fr (4.0 mm)	43 in. (109 cm)	10 ea/ca



Z8884710826E

Kangaroo™ Feeding Tubes with Standard Tip

- Weighted feeding tubes with and without flow-through stylets
- Radiopaque polyurethane material
- Centimeter markings
- MRI conditional once stylet is removed
- DEHP-free
- ENFit™ connection

Non-ENFit™ code	ENFit™ code	Includes stylet	Weight	French size	Length	Packaging
Z8884720817	Z8884720817E	No	3 g	8 Fr (2.7 mm)	36 in. (91 cm)	10 ea/ca
Z8884720841	Z8884720841E	No	3 g	8 Fr (2.7 mm)	43 in. (109 cm)	10 ea/ca
Z8884720825	Z8884720825E	Yes	3 g	8 Fr (2.7 mm)	36 in. (91 cm)	10 ea/ca
Z8884720858	Z8884720858E	Yes	3 g	8 Fr (2.7 mm)	43 in. (109 cm)	10 ea/ca
Z8884720855	Z8884720855E	Yes	3 g	8 Fr (2.7 mm)	55 in. (140 cm)	10 ea/ca
Z8884710867	Z8884710867E	Yes	5 g	8 Fr (2.7 mm)	43 in. (109 cm)	10 ea/ca
Z8884721088	Z8884721088E	Yes	5 g	10 Fr (3.3 mm)	43 in. (109 cm)	10 ea/ca
Z8884721096	Z8884721096E	No	5 g	10 Fr (3.3 mm)	43 in. (109 cm)	10 ea/ca
Z8884721055	Z8884721055E	Yes	5 g	10 Fr (3.3 mm)	55 in. (140 cm)	10 ea/ca
Z8884721211	Z8884721211E	No	5 g	12 Fr (4.0 mm)	36 in. (91 cm)	10 ea/ca
Z8884721237	Z8884721237E	No	5 g	12 Fr (4.0 mm)	43 in. (109 cm)	10 ea/ca
Z8884721252	Z8884721252E	Yes	5 g	12 Fr (4.0 mm)	43 in. (109 cm)	10 ea/ca
Z8884721255	Z8884721255E	Yes	5 g	12 Fr (4.0 mm)	55 in. (140 cm)	10 ea/ca



Z8884721211

Kangaroo™ Non-Weighted Feeding Tubes


- Clear polyurethane material
- Proximal Y-site connector with catheter tip adapter (not available in ENFit™)
- Large exit holes and open distal end
- Radiopaque stripe
- DEHP-free

Item No.	French size	Length	Packaging
Z8884711568	8 Fr (2.7 mm)	36 in. (91 cm)	10 ea/ca
Z8884711550	10 Fr (3.3 mm)	36 in. (91 cm)	10 ea/ca
Z8884711501	12 Fr (4.0 mm)	36 in. (91 cm)	10 ea/ca
Z8884711519	14 Fr (4.7 mm)	36 in. (91 cm)	10 ea/ca

Specialty Tubes

Kangaroo™ Naso-Jejunal Feeding Tubes

- Designed for precision placement in the small bowel using either endoscopy and/or fluoroscopy (includes guidewire and transfer cannula)
- Unique hollow-bolus weighted tip
- Radiopaque polyurethane material
- MRI conditional once stylet is removed
- ENFit™ connection




Z8884752505

Non-ENFit™ code	ENFit™ code	Weight	French size	Length	Packaging
Z8884752505	Z8884752505E	7 g	12 Fr (4.0 mm)	60 in. (152 cm)	2 ea/ca

Kangaroo™ Naso-Jejunal Feeding and Gastric Decompression Tubes

- For “push” insertion over guidewire (kit included)
- Radiopaque polyurethane material
- Unique adjustable-length feeding tube design
- Spiral feeding outlet holes
- Jejunal tube can be retracted within gastric portion
- MRI conditional once stylet is removed
- ENFit™ connection



Z8884752513

Non-ENFit™ code	ENFit™ code	Overall diameter	Jejunal tube diameter	Overall length	Packaging
Z8884752513	Z8884752513E	16 Fr (5.3 mm)	9 Fr (3.0 mm)	67 in. (170 cm)	2 ea/ca

Neonatal and Pediatric Feeding Tubes

Kangaroo™ and Argyle™ neonatal and pediatric feeding tubes with ENFit™ connection help address the 2014 JCAHO Tubing Misconnections Sentinel Event Alert. The ENFit™ feeding tubes:

- Are compatible with ENFit™ syringes and extension sets
- Are DEHP-free and latex-free
- Feature markings from 5-25 cm, which provide a visual indicator to assist in placement and migration checks
- Are transparent with Sentinel Line™ and Sentinel Eye™ radiopaque markings to assist in proper placement



Z461412E

Kangaroo™ Polyurethane Feeding Tubes

Non-ENFit™ code	ENFit™ code	French size	Length	Packaging
Z461355	Z461355E	3.5 Fr (1.2 mm)	12 in. (31 cm)	10 ea/ca
Z461412	Z461412E	5 Fr (1.7 mm)	20 in. (51 cm)	10 ea/ca
Z461503	Z461503E	5 Fr (1.7 mm)	36 in (91 cm)	10 ea/ca
Z461420	Z461420E	6.5 Fr (2.2 mm)	20 in. (51 cm)	10 ea/ca
Z461438	Z461438E	6.5 Fr (2.2 mm)	36 in. (91 cm)	10 ea/ca
Z461701	Z461701E	8 Fr (2.7 mm)	20 in. (51 cm)	10 ea/ca
Z461800	Z461800E	8 Fr (2.7 mm)	42 in. (107 cm)	10 ea/ca
Z461042	Z461042E	10 Fr (3.3 mm)	42 in. (107 cm)	10 ea/ca

Argyle™ PVC Feeding Tubes

Non-ENFit™ code	ENFit™ code	French size	Length	Packaging
Z461206	Z461206E	3.5 Fr (1.2 mm)	12 in. (31 cm)	10 ea/ca
Z460802	Z460802E	5 Fr (1.7 mm)	16 in. (41 cm)	50 ea/ca
Z461008	Z461008E	5 Fr (1.7 mm)	36 in. (91 cm)	50 ea/ca
Z460711	Z460711E	6.5 Fr (2.2 mm)	16 in. (41 cm)	50 ea/ca
Z460604	Z460604E	8 Fr (2.7 mm)	16 in. (41 cm)	50 ea/ca
Z460406	Z460406E	8 Fr (2.7 mm)	42 in. (107 cm)	50 ea/ca
Z460208	Z460208E	10 Fr (3.3 mm)	42 in. (107 cm)	50 ea/ca

Neonatal and Pediatric Feeding System

Kangaroo™ Feeding Tube with ENFit™ Connection Extension Sets

- PVC tubing
- DEHP-free
- With or without bifurcation (second port for irrigation and medicine delivery)
- Latex-free



Z60ENY

Non-ENFit™ code	ENFit™ code	Description	Length	Packaging
Z35ES	Z35ENS	Extension Set	35 in. (89 cm)	50 ea/ca
Z60ES	Z60ENS	Extension Set	60 in. (152 cm)	50 ea/ca
Z35EY	Z35ENY	Bifurcated Extension Set	35 in. (89 cm)	50 ea/ca
Z60EY	Z60ENY	Bifurcated Extension Set	60 in. (152 cm)	50 ea/ca

Kangaroo™ Gravity Feeding System with ENFit™ Connection

Item No.	Description	Length	Packaging
Z4VENS	Gravity Extension Set	4 in. (10.2 cm)	10 ea/ca
Z60VENS	Gravity Extension Set	60 in. (152 cm)	10 ea/ca
Z60VENY	Gravity Extension Set, Bifurcated	60 in. (152 cm)	10 ea/ca

Kangaroo™ Milk Straw with ENFit™ Connection

Item No.	Description	Length	Packaging
Z461000E	Kangaroo™ Milk Straw	5 in. (12.7 cm)	100 ea/ca



Gastrostomy Tubes

Kangaroo™ Gastrostomy Feeding Tubes with Y-Ports



- Medical-grade silicone tubing
- Rounded skin disc for external retention with feet
- Graduated shaft
- Gastric balloon for internal retention
- Rounded flush tip with open distal end
- DEHP-free

Non-ENFit™ code	ENFit™ code	French size	Balloon size	Packaging
Z8884715122	Z8884715122E	12 Fr (4.0 mm)	5 cc	5 ea/ca
Z8884715148	Z8884715148E	14 Fr (4.7 mm)	5 cc	5 ea/ca
Z8884720163	Z8884720163E	16 Fr (5.3 mm)	20 cc	5 ea/ca
Z8884720189	Z8884720189E	18 Fr (6.0 mm)	20 cc	5 ea/ca
Z8884720205	Z8884720205E	20 Fr (6.7 mm)	20 cc	5 ea/ca
Z8884720221	Z8884720221E	22 Fr (7.3 mm)	20 cc	5 ea/ca
Z8884720247	Z8884720247E	24 Fr (8.0 mm)	20 cc	5 ea/ca
Z8884720265	Z8884720265E	26 Fr (8.7 mm)	20 cc	5 ea/ca
Z8884720285	Z8884720285E	28 Fr (9.3 mm)	20 cc	5 ea/ca

Kangaroo™ Y-Port-G-Tube/P.E.G. Adapter

- For use with the Kangaroo™ P.E.G. kits

Non-ENFit™ code	ENFit™ code	Description	Packaging
Z8884751622	Z8884751622E	16 Fr Y-port/P.E.G. Adapter, Non-sterile	20 ea/ca
Z8884752042	Z8884752042E	20 Fr Y-port/P.E.G. Adapter, Non-sterile	20 ea/ca

Safety P.E.G. Kits

Kangaroo™ Safety P.E.G. Kits

- Radiopaque polyurethane material
- Non-endoscopic removal
- Easy to read placement marks
- Hydromer™* hydrophilic grafting on the “Pull” dilator tip
- DEHP-free
- Sterile

Item No.	French size	Description	Packaging
Safety Percutaneous Endoscopic Gastrostomy “Pull” Kits			
Z8884751630	16 Fr (5.3 mm)	Kit Contains: 21 G x 1-1/2 in. Needle, 25 G x 5/8 in. Needle, Dual Port Enteral Connector, 6 cc Syringe, #11 Scalpel, Catheter Assembly, 18 G Introducer Needle, External Retention Band, Looping Wire/ Guidewire, Fenestrated Drape, External Retention Disc.	2 ea/ca
Z8884752067	20 Fr (6.7 mm)	Kit Contains: 21G x 1-1/2 in. Needle, 25 G x 5/8 in. Needle, Dual Port Enteral Connector, 6 cc Syringe, #11 Scalpel, Catheter Assembly, 18 G Introducer Needle, External Retention Band, Looping Wire/ Guidewire, Fenestrated Drape, External Retention Disc.	2 ea/ca

Skin Level Tubes



Z712080

Kangaroo™ Skin Level Balloon Gastrostomy Kits

- 100% silicone material
- Durable balloon construction
- Raised feet design
- ClipStar™ enteral security clip
- DEHP-free

Contents of Kit:

Skin Level Gastrostomy Device, Step-Up Shim, Bolus Feeding Set, Continuous Feeding Set with ClipStar™ Enteral Security Clip, 60 cc Irrigation Syringe, 12 cc Syringe, Gauze Pads

Item No.	Description	Packaging	Item No.	Description	Packaging
Z712080	Balloon SLG 12 Fr x 0.8 cm	1 ea/ca	Z714350	Balloon SLG 14 Fr x 3.5 cm	1 ea/ca
Z712100	Balloon SLG 12 Fr x 1.0 cm	1 ea/ca	Z714400	Balloon SLG 14 Fr x 4.0 cm	1 ea/ca
Z712120	Balloon SLG 12 Fr x 1.2 cm	1 ea/ca	Z714450	Balloon SLG 14 Fr x 4.5 cm	1 ea/ca
Z712150	Balloon SLG 12 Fr x 1.5 cm	1 ea/ca	Z714500	Balloon SLG 14 Fr x 5.0 cm	1 ea/ca
Z712170	Balloon SLG 12 Fr x 1.7 cm	1 ea/ca	Z716080	Balloon SLG 16 Fr x 0.8 cm	1 ea/ca
Z712200	Balloon SLG 12 Fr x 2.0 cm	1 ea/ca	Z716100	Balloon SLG 16 Fr x 1.0 cm	1 ea/ca
Z712230	Balloon SLG 12 Fr x 2.3 cm	1 ea/ca	Z716120	Balloon SLG 16 Fr x 1.2 cm	1 ea/ca
Z712250	Balloon SLG 12 Fr x 2.5 cm	1 ea/ca	Z716150	Balloon SLG 16 Fr x 1.5 cm	1 ea/ca
Z712270	Balloon SLG 12 Fr x 2.7 cm	1 ea/ca	Z716170	Balloon SLG 16 Fr x 1.7 cm	1 ea/ca
Z712300	Balloon SLG 12 Fr x 3.0 cm	1 ea/ca	Z716200	Balloon SLG 16 Fr x 2.0 cm	1 ea/ca
Z712350	Balloon SLG 12 Fr x 3.5 cm	1 ea/ca	Z716230	Balloon SLG 16 Fr x 2.3 cm	1 ea/ca
Z712400	Balloon SLG 12 Fr x 4.0 cm	1 ea/ca	Z716250	Balloon SLG 16 Fr x 2.5 cm	1 ea/ca
Z712450	Balloon SLG 12 Fr x 4.5 cm	1 ea/ca	Z716270	Balloon SLG 16 Fr x 2.7 cm	1 ea/ca
Z712500	Balloon SLG 12 Fr x 5.0 cm	1 ea/ca	Z716300	Balloon SLG 16 Fr x 3.0 cm	1 ea/ca
Z714080	Balloon SLG 14 Fr x 0.8 cm	1 ea/ca	Z716350	Balloon SLG 16 Fr x 3.5 cm	1 ea/ca
Z714100	Balloon SLG 14 Fr x 1.0 cm	1 ea/ca	Z716400	Balloon SLG 16 Fr x 4.0 cm	1 ea/ca
Z714120	Balloon SLG 14 Fr x 1.2 cm	1 ea/ca	Z716450	Balloon SLG 16 Fr x 5.0 cm	1 ea/ca
Z714150	Balloon SLG 14 Fr x 1.5 cm	1 ea/ca	Z716500	Balloon SLG 16 Fr x 5.0 cm	1 ea/ca
Z714170	Balloon SLG 14 Fr x 1.7 cm	1 ea/ca	Z718080	Balloon SLG 18 Fr x 0.8 cm	1 ea/ca
Z714200	Balloon SLG 14 Fr x 2.0 cm	1 ea/ca	Z718100	Balloon SLG 18 Fr x 1.0 cm	1 ea/ca
Z714230	Balloon SLG 14 Fr x 2.3 cm	1 ea/ca	Z718120	Balloon SLG 18 Fr x 1.2 cm	1 ea/ca
Z714250	Balloon SLG 14 Fr x 2.5 cm	1 ea/ca	Z718150	Balloon SLG 18 Fr x 1.5 cm	1 ea/ca
Z714270	Balloon SLG 14 Fr x 2.7 cm	1 ea/ca	Z718170	Balloon SLG 18 Fr x 1.7 cm	1 ea/ca
Z714300	Balloon SLG 14 Fr x 3.0 cm	1 ea/ca	Z718200	Balloon SLG 18 Fr x 2.0 cm	1 ea/ca

Kangaroo™ Skin Level Balloon Gastrostomy Kits (continued)

Item No.	Description	Packaging	Item No.	Description	Packaging
Z718230	Balloon SLG 18 Fr x 2.3 cm	1 ea/ca	Z720350	Balloon SLG 20 Fr x 3.5 cm	1 ea/ca
Z718250	Balloon SLG 18 Fr x 2.5 cm	1 ea/ca	Z720400	Balloon SLG 20 Fr x 4.0 cm	1 ea/ca
Z718270	Balloon SLG 18 Fr x 2.7 cm	1 ea/ca	Z720450	Balloon SLG 20 Fr x 4.5 cm	1 ea/ca
Z718300	Balloon SLG 18 Fr x 3.0 cm	1 ea/ca	Z720500	Balloon SLG 20 Fr x 5.0 cm	1 ea/ca
Z718350	Balloon SLG 18 Fr x 3.5 cm	1 ea/ca	Z724080	Balloon SLG 24 Fr x 0.8 cm	1 ea/ca
Z718400	Balloon SLG 18 Fr x 4.0 cm	1 ea/ca	Z724100	Balloon SLG 24 Fr x 1.0 cm	1 ea/ca
Z718450	Balloon SLG 18 Fr x 4.5 cm	1 ea/ca	Z724120	Balloon SLG 24 Fr x 1.2 cm	1 ea/ca
Z718500	Balloon SLG 18 Fr x 5.0 cm	1 ea/ca	Z724150	Balloon SLG 24 Fr x 1.5 cm	1 ea/ca
Z720080	Balloon SLG 20 Fr x 0.8 cm	1 ea/ca	Z724170	Balloon SLG 24 Fr x 1.7 cm	1 ea/ca
Z720100	Balloon SLG 20 Fr x 1.0 cm	1 ea/ca	Z724200	Balloon SLG 24 Fr x 2.0 cm	1 ea/ca
Z720120	Balloon SLG 20 Fr x 1.2 cm	1 ea/ca	Z724230	Balloon SLG 24 Fr x 2.3 cm	1 ea/ca
Z720150	Balloon SLG 20 Fr x 1.5 cm	1 ea/ca	Z724250	Balloon SLG 24 Fr x 2.5 cm	1 ea/ca
Z720170	Balloon SLG 20 Fr x 1.7 cm	1 ea/ca	Z724270	Balloon SLG 24 Fr x 2.7 cm	1 ea/ca
Z720200	Balloon SLG 20 Fr x 2.0 cm	1 ea/ca	Z724300	Balloon SLG 24 Fr x 3.0 cm	1 ea/ca
Z720230	Balloon SLG 20 Fr x 2.3 cm	1 ea/ca	Z724350	Balloon SLG 24 Fr x 3.5 cm	1 ea/ca
Z720250	Balloon SLG 20 Fr x 2.5 cm	1 ea/ca	Z724400	Balloon SLG 24 Fr x 4.0 cm	1 ea/ca
Z720270	Balloon SLG 20 Fr x 2.7 cm	1 ea/ca	Z724450	Balloon SLG 24 Fr x 4.5 cm	1 ea/ca
Z720300	Balloon SLG 20 Fr x 3.0 cm	1 ea/ca	Z724500	Balloon SLG 24 Fr x 5.0 cm	1 ea/ca

Accessories

Item No.	Description	Packaging
Z8884741805	Stoma Measuring Device, Sterile	5 ea/ca
Z8884741813	Bolus Feeding Set, Non-sterile	10 ea/ca
Z8884741821	Continuous Feeding Set, Non-sterile	10 ea/ca



Standard Syringes

Monoject™ Oral Syringes – Standard Tip

Item No.	Description	Packaging
Z401S	1 mL Non-sterile Purple Oral Syringe	100 ea/ca
Z403S	3 mL Non-sterile Purple Oral Syringe	500 ea/ca
Z406S	6 mL Non-sterile Purple Oral Syringe	500 ea/ca
Z410S	10 mL Non-sterile Purple Oral Syringe	500 ea/ca
Z401SG	1 mL Sterile Purple Oral Syringe	240 ea/ca
Z403SG	3 mL Sterile Purple Oral Syringe	480 ea/ca
Z406SG	6 mL Sterile Purple Oral Syringe	400 ea/ca
Z412SG	12 mL Sterile Purple Oral Syringe	192 ea/ca
Z435SG	35 mL Sterile Purple Oral Syringe	160 ea/ca
Z460SG	60 mL Sterile Purple Oral Syringe	120 ea/ca



ENFit™ Syringes

Monoject™ Enteral Syringes – ENFit™ Connection

Item No.	Description	Packaging
Z8881101015	1 mL Non-sterile Enteral Syringe	3200 ea/ca
Z8881103015	3 mL Non-sterile Enteral Syringe	2000 ea/ca
Z8881106015	6 mL Non-sterile Enteral Syringe	1200 ea/ca
Z8881112015	12 mL Non-sterile Enteral Syringe	1000 ea/ca
Z8881135015	35 mL Non-sterile Enteral Syringe	300 ea/ca
Z8881160015	60 mL Non-sterile Enteral Syringe	175 ea/ca
Z401SE	1 mL Sterile Enteral Syringe	240 ea/ca
Z403SE	3 mL Sterile Enteral Syringe	400 ea/ca
Z406SE	6 mL Sterile Enteral Syringe	400 ea/ca
Z412SE	12 mL Sterile Enteral Syringe	480 ea/ca
Z435SE	35 mL Sterile Enteral Syringe	160 ea/ca
Z460SE	60 mL Sterile Enteral Syringe	120 ea/ca

Monoject™ Enteral Syringe Tip Caps with ENFit™ Connection

Item No.	Description	Packaging
Z400SE	ENFit™ Syringe Tip Cap, Sterile	1000 ea/ca
Z400BE	ENFit™ Syringe Tip Cap, Non-sterile	1000 ea/ca




Z400BE

Monoject™ Medicine Straw with ENFit™ Connection

Item No.	Description	Length	Packaging
Z462000E	Monoject™ Medicine Straw	2 in (5.08 cm)	100 ea/ca
Z464000E	Monoject™ Medicine Straw	4 in (10.16 cm)	100 ea/ca

Monoject™ UV/Light Inhibiting Transfer Bag

Item No.	Description	Packaging
Z46312	UV/Light Inhibiting Transfer Bag	250 ea/ca



Z46312

Monoject™ Universal Bottle Adapter with ENFit™ Connection


Item No.	Description	Packaging
Z4699E	Universal Bottle Adapter	50 ea/ca



Z4699E

Salem Sump™ Tubes

Salem Sump™ PVC Tubes




Z8888268060

- Double lumen PVC tube—one for suction drainage, one for sump vent
- Centimeter markings assist the clinician with tube placement by providing reference points on the tube and allow for tracking of tube migration after tube securement
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Funnel end of vent lumen can be used to cap 5-in-1 adapter
- Sterile

Item No.	French size	Length	Packaging
Z8888268060	6 Fr (2.0 mm)	24 in. (61 cm)	10 ea/ca
Z8888268086	8 Fr (2.7 mm)	24 in. (61 cm)	10 ea/ca
Z8888264911	10 Fr (3.3 mm)	36 in. (91 cm)	50 ea/ca
Z8888264929	12 Fr (4.0 mm)	48 in. (122 cm)	50 ea/ca
Z8888264945	14 Fr (4.7 mm)	48 in. (122 cm)	50 ea/ca
Z8888264960	16 Fr (5.3 mm)	48 in. (122 cm)	50 ea/ca
Z8888264986	18 Fr (6.0 mm)	48 in. (122 cm)	50 ea/ca

Salem Sump™ PVC Tubes with Anti-Reflux Valve



Z8888266106

- Double lumen PVC tube—one for suction drainage, one for sump vent
- Gastric reflux prevention valve
- Centimeter markings assist the clinician with tube placement by providing reference points on the tube and allow for tracking of tube migration after tube securement
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Funnel end of vent lumen can be used to cap 5-in-1 adapter
- Sterile

Item No.	French size	Length	Packaging
Z8888266106	10 Fr (3.3 mm)	36 in. (91 cm)	50 ea/ca
Z8888266114	12 Fr (4.0 mm)	48 in. (122 cm)	50 ea/ca
Z8888266122	14 Fr (4.7 mm)	48 in. (122 cm)	50 ea/ca
Z8888266130	16 Fr (5.3 mm)	48 in. (122 cm)	50 ea/ca
Z8888266148	18 Fr (6.0 mm)	48 in. (122 cm)	50 ea/ca

Salem Sump™ Anti-Reflux Valves

- Gastric reflux prevention valve
- Sterile
- Allows for maximum air flow, assists to reduce stomach mucosa damage

Item No.	Description	Packaging
Z8888266197	Valve only	10 ea/ca




Z8888266197

Salem Sump™ Silicone Tubes

- Double lumen silicone tube—one for suction drainage, one for venting action
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Removable 5-in-1 connector
- Capping system with sump vent lumen cap, suction lumen adapter and suction lumen adapter cap
- Sterile

Item No.	French size	Length	Packaging
Z8888265108	10 Fr (3.3 mm)	36 in. (91 cm)	10 ea/ca
Z8888265116	12 Fr (4.0 mm)	48 in. (122 cm)	10 ea/ca
Z8888265124	14 Fr (4.7 mm)	48 in. (122 cm)	10 ea/ca
Z8888265132	16 Fr (5.3 mm)	48 in. (122 cm)	10 ea/ca
Z8888265140	18 Fr (6.0 mm)	48 in. (122 cm)	10 ea/ca




Z8888265108

Salem Sump™ Silicone Tubes with Anti-Reflux Valve

- Double lumen silicone tube—one for suction drainage, one for venting action
- Gastric reflux prevention valve
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Removable 5-in-1 connector
- Capping system with sump vent lumen cap, suction lumen adapter and suction lumen adapter cap
- Sterile


Item No.	French size	Length	Packaging
Z8888266213	10 Fr (3.3 mm)	36 in. (91 cm)	10 ea/ca
Z8888266221	12 Fr (4.0 mm)	48 in. (122 cm)	10 ea/ca
Z8888266247	14 Fr (4.7 mm)	48 in. (122 cm)	10 ea/ca
Z8888266262	16 Fr (5.3 mm)	48 in. (122 cm)	10 ea/ca
Z8888266288	18 Fr (6.0 mm)	48 in. (122 cm)	10 ea/ca



Z8888266213

Stomach Tubes

Argyle™ Stomach Tubes (Levin Type)



Z8888265314

- Clear, thermosensitive PVC tube with four eyes
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Open distal end
- Integral funnel connector
- Markings 18 in. (46 cm), 22 in. (56 cm), 26 in. (66 cm) and 30 in. (76 cm) from distal end
- Non-sterile

Item No.	French size	Length	Packaging
Z8888265314	12 Fr (4.0 mm)	48 in. (122 cm)	50 ea/ca
Z8888265322	16 Fr (5.3 mm)	48 in. (122 cm)	50 ea/ca
Z8888265330	18 Fr (6.0 mm)	48 in. (122 cm)	50 ea/ca



Z155712

Stomach Tubes (Levin Type)

- Tube has four staggered eyes and center opening to reduce trauma
- Depth marks begin 20 in. from distal tip with three additional marks at 4 in. intervals
- Smooth, rounded tip reduces insertion trauma
- Non-sterile

Item No.	French size	Length	Packaging
Z155710	12 Fr (4.0 mm)	48 in. (122 cm)	50 ea/ca
Z155711	14 Fr (4.7 mm)	48 in. (122 cm)	50 ea/ca
Z155712	16 Fr (5.3 mm)	48 in. (122 cm)	50 ea/ca
Z155713	18 Fr (6.0 mm)	48 in. (122 cm)	50 ea/ca

Argyle™ Edlich Gastric Lavage Trays


- Large 140 cc (5 oz) graduated Monoject™ syringes
- Transparent PVC gastric tube with closed end and four eyes
- Depth marks at 18 in. (46 cm), 22 in. (56 cm), 26 in. (66 cm) and 30 in. (76 cm)
- Clamp

Item No.	French size	Length	Packaging
Z8888750018	Large, 34 Fr (11.3 mm)	36 in. (91 cm)	20 ea/ca
Z8888750513	Medium/Pediatric, 24 Fr (8.0 mm)	36 in. (91 cm)	20 ea/ca



Z8888750018





Z7771810E

Kangaroo™ Multifunctional Ports with Safe Enteral Connections

- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- All-in-one closed port system
- Allows for feeding, suction, irrigation and medication
- Internal ARV design
- DEHP-free
- Radiopaque tubing

Non-ENFit™ code	ENFit™ code	Description	Packaging
Z771111	Z771111E	Kangaroo™ Multifunctional Port with Single Lumen Adapter	10 ea/ca

Salem Sump™ Tubes

Non-ENFit™ code	ENFit™ code	Description	Packaging
Z7771010	Z7771010E	10 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10 ea/ca
Z7771210	Z7771210E	12 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10 ea/ca
Z7771410	Z7771410E	14 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10 ea/ca
Z7771610	Z7771610E	16 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10 ea/ca
Z7771810	Z7771810E	18 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	10 ea/ca

CO₂ Detectors

Kangaroo™ CO₂ Detector



Z777702

- Proprietary colorimetric paper technology
- Detects presence of CO₂ within seconds
- Easily attaches to either a bellows or syringe
- Requires no change in insertion technique
- Convenient bellows device:
 - Clears obstructions from distal tip
 - Prevents CO₂ from entering tube prematurely
 - Facilitates airflow for quick results

Non-ENFit™ code	ENFit™ code	Description	Packaging
Z777702	Z777702E	Kangaroo™ CO ₂ Detector	10 ea/ca



Index

Code	Product	Page	Code	Product	Page
Z1KCP	Kangaroo™ Connect ENPlus Spike ENtelliSet with ENFit™ Connector, Non-sterile	7	Z412SE	12 mL Sterile Enteral Syringe	20
Z1000KCP	Kangaroo™ Connect ENTelliset with ENFit™ Connector, 1000 mL, Non-sterile	7	Z435SG	35 mL Sterile Purple Oral Syringe	20
Z1054715	Kangaroo™ ePump Battery	4	Z435SE	35 mL Sterile Enteral Syringe	20
Z155710	12 Fr, 48 in. (122 cm) Stomach Tubes, Levin Type	24	Z460SE	60 mL Sterile Enteral Syringe	20
Z155711	14 Fr, 48 in. (122 cm) Stomach Tubes, Levin Type	24	Z460SG	60 mL Sterile Purple Oral Syringe	20
Z155712	16 Fr, 48 in. (122 cm) Stomach Tubes, Levin Type	24	Z460208/ Z460208E	10 Fr, 42 in. (107 cm) Argyle™ PVC Feeding Tube	14
Z155713	18 Fr, 48 in. (122 cm) Stomach Tubes, Levin Type	24	Z460406/ Z460406E	8 Fr, 42 in. (107 cm) Argyle™ PVC Feeding Tube	14
Z35ES/Z35ENS	35 in. Extension Set	14	Z460604/ Z460604E	8 Fr, 16 in. (41 cm) Argyle™ PVC Feeding Tube	14
Z35EY/Z35ENY	35 in. Bifurcated Extension Set	14	Z460711/ Z460711E	6.5 Fr, 16 in. (41 cm) Argyle™ PVC Feeding Tube	14
Z382491	Kangaroo™ ePump Power Adapter	4	Z460802/ Z460802E	5 Fr, 16 in. (41 cm) Argyle™ PVC Feeding Tube	14
Z382492	Kangaroo™ ePump Pole Clamp	4	Z460836E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 36 in. (91 cm)	10
Z383400	Kangaroo™ Joey Enteral Feeding Pump with Pole Clamp	5	Z460843E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 43 in. (109 cm)	10
Z383491	Kangaroo™ Joey Power Cord	5	Z460855E	Kangaroo™ Feeding Tube with IRIS Technology 8 Fr x 55 in. (140 cm)	10
Z383493	Kangaroo™ Joey Pole Clamp	5	Z461000E	Kangaroo™ Milk Straw	15
Z384400	Kangaroo™ Connect Enteral Feeding Pump with Pole Clamp and Power Adapter	6	Z461008/ Z461008E	5 Fr, 36 in. (91 cm) Argyle™ PVC Feeding Tube	14
Z384491	Kangaroo™ Connect Power Cord with Adapter	6	Z461036E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 36 in. (91 cm)	10
Z384492	Kangaroo™ Connect Pole Clamp	6	Z461043E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 43 in. (109 cm)	10
Z384494	Kangaroo™ Connect Car Charger	6	Z461055E	Kangaroo™ Feeding Tube with IRIS Technology 10 Fr x 55 in. (140 cm)	10
Z384500CBW	Kangaroo™ Connect WiFi/Ethernet Communications Hub	7	Z461236E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 36 in (91 cm)	10
Z384500CBC	Kangaroo™ Connect Cellular Communications Hub	7	Z461243E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 43 in. (109 cm)	10
Z386100	Kangaroo™ IRIS Console with Accessories (Global Power Adapter, Battery, Interface Cable and Pole Clamp)	10	Z461255E	Kangaroo™ Feeding Tube with IRIS Technology 12 Fr x 55 in. (140 cm)	10
Z388101	Kangaroo™ IRIS Power Adapter	10	Z461042/ Z461042E	10 Fr, 42 in. (107 cm) Kangaroo™ Polyurethane Feeding Tube	13
Z388102	Kangaroo™ IRIS Replacement Battery	10	Z461206/ Z461206E	3.5 Fr, 12 in. (31 cm) Argyle™ PVC Feeding Tube	14
Z388103	Kangaroo™ IRIS Interface Cable	10	Z461355/ Z461355E	3.5 Fr, 12 in. (31 cm) Kangaroo™ Polyurethane Feeding Tube	13
Z388104	Kangaroo™ IRIS Mounting Clamp	10	Z461412/ Z461412E	5 Fr, 20 in. (51 cm) Kangaroo™ Polyurethane Feeding Tube	13
Z388105E	Kangaroo™ IRIS Insufflation Device	10	Z461420/ Z461420E	6.5 Fr, 20 in. (51 cm) Kangaroo™ Polyurethane Feeding Tube	13
Z4VENS	4 in. Gravity Extension Set	14	Z461438/ Z461438E	6.5 Fr, 36 in. (91 cm) Kangaroo™ Polyurethane Feeding Tube	13
Z400SE	ENFit™ Syringe Tip Cap, Sterile	21	Z461503/ Z461503E	5 Fr, 36 in. (91 cm) Kangaroo™ Polyurethane Feeding Tube	13
Z400BE	ENFit™ Syringe Tip Cap, Non-sterile	21	Z461701/ Z461701E	8 Fr, 20 in. (51 cm) Kangaroo™ Polyurethane Feeding Tube	13
Z401S	1 mL Non-sterile Purple Oral Syringe	20			
Z401SE	1 mL Sterile Enteral Syringe	20			
Z401SG	1 mL Sterile Purple Oral Syringe	20			
Z403S	3 mL Non-sterile Purple Oral Syringe	20			
Z403SE	3 mL Sterile Enteral Syringe	20			
Z403SG	3 mL Sterile Purple Oral Syringe	20			
Z406S	6 mL Non-sterile Purple Oral Syringe	20			
Z406SE	6 mL Sterile Enteral Syringe	20			
Z406SG	6 mL Non-sterile Purple Oral Syringe	20			
Z410S	10 mL Non-sterile Purple Oral Syringe	20			
Z412SG	12 mL Sterile Purple Oral Syringe	20			

Code	Product	Page	Code	Product	Page
Z461800/ Z461800E	8 Fr, 42 in. (107 cm) Kangaroo™ Polyurethane Feeding Tube	13	Z716150	Balloon SLG 16 Fr x 1.5 cm	18
Z462000E	Monoject™ Medicine Straw 6 cm	21	Z716154	Kangaroo™ ePump Sterile 100 mL Burette Set	4
Z46312	UV/Light Inhibiting Transfer Bag	21	Z716170	Balloon SLG 16 Fr x 1.7 cm	18
Z464000E	Monoject™ Medicine Straw 10 cm	21	Z716200	Balloon SLG 16 Fr x 2.0 cm	18
Z4699E	Universal Bottle Adapter	21	Z716230	Balloon SLG 16 Fr x 2.3 cm	18
Z482400J	Kangaroo™ ePump with Pole Clamp	4	Z716250	Balloon SLG 16 Fr x 2.5 cm	18
Z500KCP	Kangaroo™ Connect ENTelliset with ENFit™ Connector, 500 mL, Non-sterile	7	Z716270	Balloon SLG 16 Fr x 2.7 cm	18
Z60ES/Z60ENS	60 in. Extension Set	14	Z716300	Balloon SLG 16 Fr x 3.0 cm	18
Z60EY/Z60ENY	60 in. Bifurcated Extension Set	14	Z716350	Balloon SLG 16 Fr x 3.5 cm	18
Z60VENS	60 in. Gravity Extension Set	14	Z716400	Balloon SLG 16 Fr x 4.0 cm	18
Z60VENY	60 in. Gravity Extension Set, Bifurcated	14	Z716450	Balloon SLG 16 Fr x 4.5 cm	18
Z674668	Kangaroo™ ePump 500 mL Feed Set with 500 mL Flush Bag	4	Z716500	Balloon SLG 16 Fr x 5.0 cm	18
Z702505	1000 mL Large Bore Gravity Set	9	Z718080	Balloon SLG 18 Fr x 0.8 cm	18
Z712080	Balloon SLG 12 Fr x 0.8 cm	18	Z718100	Balloon SLG 18 Fr x 1.0 cm	18
Z712100	Balloon SLG 12 Fr x 1.0 cm	18	Z718120	Balloon SLG 18 Fr x 1.2 cm	18
Z712120	Balloon SLG 12 Fr x 1.2 cm	18	Z718150	Balloon SLG 18 Fr x 1.5 cm	18
Z712150	Balloon SLG 12 Fr x 1.5 cm	18	Z718170	Balloon SLG 18 Fr x 1.7 cm	18
Z712170	Balloon SLG 12 Fr x 1.7 cm	18	Z718200	Balloon SLG 18 Fr x 2.0 cm	18
Z712200	Balloon SLG 12 Fr x 2.0 cm	18	Z718230	Balloon SLG 18 Fr x 2.3 cm	19
Z712230	Balloon SLG 12 Fr x 2.3 cm	18	Z718250	Balloon SLG 18 Fr x 2.5 cm	19
Z712250	Balloon SLG 12 Fr x 2.5 cm	18	Z718270	Balloon SLG 18 Fr x 2.7 cm	19
Z712270	Balloon SLG 12 Fr x 2.7 cm	18	Z718300	Balloon SLG 18 Fr x 3.0 cm	19
Z712300	Balloon SLG 12 Fr x 3.0 cm	18	Z718350	Balloon SLG 18 Fr x 3.5 cm	19
Z712350	Balloon SLG 12 Fr x 3.5 cm	18	Z718400	Balloon SLG 18 Fr x 4.0 cm	19
Z712400	Balloon SLG 12 Fr x 4.0 cm	18	Z718450	Balloon SLG 18 Fr x 4.5 cm	19
Z712450	Balloon SLG 12 Fr x 4.5 cm	18	Z718500	Balloon SLG 18 Fr x 5.0 cm	19
Z712500	Balloon SLG 12 Fr x 5.0 cm	18	Z720080	Balloon SLG 20 Fr x 0.8 cm	19
Z714080	Balloon SLG 14 Fr x 0.8 cm	18	Z720100	Balloon SLG 20 Fr x 1.0 cm	19
Z714100	Balloon SLG 14 Fr x 1.0 cm	18	Z720120	Balloon SLG 20 Fr x 1.2 cm	19
Z714120	Balloon SLG 14 Fr x 1.2 cm	18	Z720150	Balloon SLG 20 Fr x 1.5 cm	19
Z714150	Balloon SLG 14 Fr x 1.5 cm	18	Z720170	Balloon SLG 20 Fr x 1.7 cm	19
Z714170	Balloon SLG 14 Fr x 1.7 cm	18	Z720200	Balloon SLG 20 Fr x 2.0 cm	19
Z714200	Balloon SLG 14 Fr x 2.0 cm	18	Z720230	Balloon SLG 20 Fr x 2.3 cm	19
Z714230	Balloon SLG 14 Fr x 2.3 cm	18	Z720250	Balloon SLG 20 Fr x 2.5 cm	19
Z714250	Balloon SLG 14 Fr x 2.5 cm	18	Z720270	Balloon SLG 20 Fr x 2.7 cm	19
Z714270	Balloon SLG 14 Fr x 2.7 cm	18	Z720300	Balloon SLG 20 Fr x 3.0 cm	19
Z714300	Balloon SLG 14 Fr x 3.0 cm	18	Z720350	Balloon SLG 20 Fr x 3.5 cm	19
Z714350	Balloon SLG 14 Fr x 3.5 cm	18	Z720400	Balloon SLG 20 Fr x 4.0 cm	19
Z714400	Balloon SLG 14 Fr x 4.0 cm	18	Z720450	Balloon SLG 20 Fr x 4.5 cm	19
Z714450	Balloon SLG 14 Fr x 4.5 cm	18	Z720500	Balloon SLG 20 Fr x 5.0 cm	19
Z714500	Balloon SLG 14 Fr x 5.0 cm	18	Z724080	Balloon SLG 24 Fr x 0.8 cm	19
Z716080	Balloon SLG 16 Fr x 0.8 cm	18	Z724100	Balloon SLG 24 Fr x 1.0 cm	19
Z716100	Balloon SLG 16 Fr x 1.0 cm	18	Z724120	Balloon SLG 24 Fr x 1.2 cm	19
Z716120	Balloon SLG 16 Fr x 1.2 cm	18	Z724150	Balloon SLG 24 Fr x 1.5 cm	19
			Z724170	Balloon SLG 24 Fr x 1.7 cm	19
			Z724200	Balloon SLG 24 Fr x 2.0 cm	19

Code	Product	Page	Code	Product	Page
Z724230	Balloon SLG 24 Fr x 2.3 cm	19	Z7771210/ Z7771210E	12 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	26
Z724250	Balloon SLG 24 Fr x 2.5 cm	19	Z7771410/ Z7771410E	14 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	26
Z724270	Balloon SLG 24 Fr x 2.7 cm	19	Z7771610/ Z7771610E	16 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	26
Z724300	Balloon SLG 24 Fr x 3.0 cm	19	Z7771810/ Z7771810E	18 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	26
Z724350	Balloon SLG 24 Fr x 3.5 cm	19	Z777702/ Z777702E	Kangaroo™ CO ₂ Detector	26
Z724400	Balloon SLG 24 Fr x 4.0 cm	19	Z84065/ Z84065E	I.V. Pouch with 60 cc Piston Syringe	9
Z724450	Balloon SLG 24 Fr x 4.5 cm	19	Z8881101015	1 mL Non-sterile Enteral Syringe	20
Z724500	Balloon SLG 24 Fr x 5.0 cm	19	Z8881103015	3 mL Non-sterile Enteral Syringe	20
Z762055	Kangaroo™ Joey 500 mL Pump Set	5	Z8881106015	6 mL Non-sterile Enteral Syringe	20
Z763656	Kangaroo™ Joey 1000 mL Pump Set	5	Z8881112015	12 mL Non-sterile Enteral Syringe	20
Z763662	Kangaroo™ Joey 1000 mL Pump Set with Flush Bag	5	Z8881135015	35 mL Non-sterile Enteral Syringe	20
Z765100	Kangaroo™ Joey Safety Screw Spike 1000 mL Flush Bag	5	Z8881160015	60 mL Non-sterile Enteral Syringe	20
Z765000	Kangaroo™ ENFit™ Transition Connector	9	Z8884700108/ Z8884700108E	Irrigation Kit Syringe with 500 cc Container and Basin	9
Z765559	Kangaroo™ Joey ENPlus Spike Set	5	Z8884700116/ Z8884700116E	Irrigation Kit Syringe with 500 cc Container	9
Z77000FD	Kangaroo™ Connect ENPlus Spike Set, Non-sterile	7	Z8884700207/ Z8884700207E	Pump Set 4 ft. Tubing Extension	9
Z770025	Kangaroo™ Joey Black Mini Backpack	6	Z8884702500	1000 mL Gravity Bag	9
Z770026	Kangaroo™ Joey Blue Mini Backpack	6	Z8884705008/ Z8884705008E	Y-site Extension Set	9
Z770027	Kangaroo™ Joey Green Mini Backpack	6	Z8884710826/ Z8884710826E	8 Fr, 55 in. (140 cm) Kangaroo™ Feeding Tube with Dobbhoff™ Tip	11
Z770028	Kangaroo™ Joey Blue Super Mini Backpack	6	Z8884710842/ Z8884710842E	8 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Dobbhoff™ Tip, No Stylet	11
Z770029	Kangaroo™ Joey Pink Super Mini Backpack	6	Z8884710859/ Z8884710859E	8 Fr, 43 in. (109 cm), Kangaroo™ Feeding Tube with Dobbhoff™ Tip	11
Z770031	Kangaroo™ Joey Black Super Mini Backpack	6	Z8884710867/ Z8884710867E	8 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Standard Tip	11
Z770032	Kangaroo™ Joey Green Super Mini Backpack	6	Z8884711006/ Z8884711006E	10 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Dobbhoff™ Tip	11
Z770034	Kangaroo™ Joey Pink Mini Backpack	6	Z8884711246/ Z8884711246E	12 Fr, 43 in. (109 cm), Kangaroo™ Feeding Tube with Dobbhoff™ Tip, No Stylet	11
Z770035S	Kangaroo™ Connect Backpack Purple, Small	8	Z8884711253/ Z8884711253E	12 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Dobbhoff™ Tip	11
Z770035M	Kangaroo™ Connect Backpack Purple, Medium	8	Z8884711501	12 Fr, 36 in. (91 cm) Kangaroo™ Non-weighted Feeding Tube	12
Z770035L	Kangaroo™ Connect Backpack Purple, Large	8	Z8884711519	14 Fr, 36 in. (91 cm) Kangaroo™ Non-weighted Feeding Tube	12
Z770037S	Kangaroo™ Connect Backpack Black, Small	8	Z8884711550	10 Fr, 36 in. (91 cm) Kangaroo™ Non-weighted Feeding Tube	12
Z770037M	Kangaroo™ Connect Backpack Black, Medium	8	Z8884711568	8 Fr, 36 in. (91 cm) Kangaroo™ Non-weighted Feeding Tube	12
Z770037L	Kangaroo™ Connect Backpack Black, Large	8	Z8884715122/ Z8884715122E	12 Fr, 5 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
Z77100FD	Kangaroo™ Connect 1000 mL Bag Set, Non-sterile	7	Z8884715148/ Z8884715148E	14 Fr, 5 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16
Z771111/ Z771111E	Kangaroo™ Multifunctional Port with Single Lumen Adapter	26			
Z772055	Kangaroo™ ePump 500 mL Pump Set	4			
Z773656	Kangaroo™ ePump 1000 mL Pump Set	4			
Z773662	Kangaroo™ ePump 1000 mL Feed Set with 1000 mL Flush Bag	4			
Z77500FD	Kangaroo™ Connect 500 mL Bag Set, Non-sterile	7			
Z775100	Kangaroo™ ePump ENPlus Spike with 1000 mL Flush Bag	4			
Z775659	Kangaroo™ ePump ENPlus Spike Set	4			
Z776150	Kangaroo™ ePump 100 mL Burette Recertification Set	4			
Z7771010/ Z7771010E	10 Fr Salem Sump™ Tube with Kangaroo™ Multifunctional Port	26			

Code	Product	Page	Code	Product	Page
Z8884720163/ Z8884720163E	16 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16	Z8888264911	10 Fr, 36 in. (91 cm) Salem Sump™ PVC Tube	22
Z8884720189/ Z8884720189E	18 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16	Z8888264929	12 Fr, 48 in. (122 cm) Salem Sump™ PVC Tube	22
Z8884720205/ Z8884720205E	20 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16	Z8888264945	14 Fr, 48 in. (122 cm) Salem Sump™ PVC Tube	22
Z8884720221/ Z8884720221E	22 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16	Z8888264960	16 Fr, 48 in. (122 cm) Salem Sump™ PVC Tube	22
Z8884720247/ Z8884720247E	24 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16	Z8888264986	18 Fr, 48 in. (122 cm) Salem Sump™ PVC Tube	22
Z8884720265/ Z8884720265E	26 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16	Z8888265108	10 Fr, 36 in. (91 cm) Salem Sump™ Silicone Tube	23
Z8884720285/ Z8884720285E	28 Fr, 20 cc Kangaroo™ Gastrostomy Feeding Tube with Y-Port	16	Z8888265116	12 Fr, 48 in. (122 cm) Salem Sump™ Silicone Tube	23
Z8884720817/ Z8884720817E	8 Fr, 36 in. (91 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11	Z8888265124	14 Fr, 48 in. (122 cm) Salem Sump™ Silicone Tube	23
Z8884720825/ Z8884720825E	8 Fr, 36 in. (91 cm) Kangaroo™ Feeding Tube with Standard Tip	11	Z8888265132	16 Fr, 48 in. (122 cm) Salem Sump™ Silicone Tube	23
Z8884720841/ Z8884720841E	8 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11	Z8888265140	18 Fr, 48 in. (122 cm) Salem Sump™ Silicone Tube	23
Z8884720855/ Z8884720855E	8 Fr, 55 in. (140 cm) Kangaroo™ Feeding Tube with Standard Tip	11	Z8888265314	12 Fr, 48 in. (122 cm) Argyle™ Stomach Tube, Levin Type	24
Z8884720858/ Z8884720858E	8 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Standard Tip	11	Z8888265322	16 Fr, 48 in. (122 cm) Argyle™ Stomach Tube, Levin Type	24
Z8884721055/ Z8884721055E	10 Fr, 55 in. (140 cm) Kangaroo™ Feeding Tube with Standard Tip	11	Z8888265330	18 Fr, 48 in. (122 cm) Argyle™ Stomach Tube, Levin Type	24
Z8884721088/ Z8884721088E	10 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Standard Tip	11	Z8888266106	10 Fr, 36 in. (91 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	22
Z8884721096/ Z8884721096E	10 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11	Z8888266114	12 Fr, 48 in. (122 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	22
Z8884721211/ Z8884721211E	12 Fr, 36 in. (91 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11	Z8888266122	14 Fr, 48 in. (122 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	22
Z8884721237/ Z8884721237E	12 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Standard Tip, No Stylet	11	Z8888266130	16 Fr, 48 in. (122 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	22
Z8884721252/ Z8884721252E	12 Fr, 43 in. (109 cm) Kangaroo™ Feeding Tube with Standard Tip	11	Z8888266148	18 Fr, 48 in. (122 cm) Salem Sump™ PVC Tube with Anti-Reflux Valve	22
Z8884721255/ Z8884721255E	12 Fr, 55 in. (140 cm) Kangaroo™ Feeding Tube with Standard Tip	11	Z8888266197	Salem Sump™ Anti-Reflux Valve	23
Z8884741805	Stoma Measuring Device, Sterile	19	Z8888266213	10 Fr, 36 in. (91 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	23
Z8884741813	Bolus Feeding Set, Non-sterile	19	Z8888266221	12 Fr, 48 in. (122 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	23
Z8884741821	Continuous Feeding Set, Non-sterile	19	Z8888266247	14 Fr, 48 in. (122 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	23
Z8884751622/ Z8884751622E	16 Fr Y-port/P.E.G. Adapter, Non-sterile	16	Z8888266262	16 Fr, 48 in. (122 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	23
Z8884751630	16 Fr, Safety Percutaneous Endoscopic Gastrostomy "Pull" Kits	17	Z8888266288	18 Fr, 48 in. (122 cm) Salem Sump™ Silicone Tube with Anti-Reflux Valve	23
Z8884752042/ Z8884752042E	20 Fr Y-port/P.E.G. Adapter, Non-sterile	16	Z8888268060	6 Fr, 24 in. (61 cm) Salem Sump™ PVC Tub	22
Z8884752067	20 Fr, Safety Percutaneous Endoscopic Gastrostomy "Pull" Kits	17	Z8888268086	8 Fr, 24 in. (61 cm) Salem Sump™ PVC Tube	22
Z8884752505/ Z8884752505E	12 Fr, 60 in. (152 cm) Kangaroo™ Naso-Jejunal Feeding Tube	12	Z8888750018	Argyle™ Edlich Gastric Lavage Trays, Large	25
Z8884752513/ Z8884752513E	16 Fr, 67 in. (170 cm) Kangaroo™ Naso-Jejunal Feeding and Gastric Decompression Tube	12	Z8888750513	Argyle™ Edlich Gastric Lavage Trays, Medium	25
			ZF010506	Kangaroo™ Joey Battery	5

For more information or to order, contact your local
Cardinal Health Canada **sales representative today**
or **call 1.888.291.5033**

The International Organization for Standardization (ISO) and the Global Enteral Device Supplier Association (GEDSA) have created a new standard design for enteral connectors known as ENFit™ connectors. ENFit™ connectors are intended to improve patient safety and decrease the risk of placing IV medicines in an enteral feeding device.

© 2019 Cardinal Health Canada. All Rights Reserved. CARDINAL HEALTH, Kangaroo, Salem Sump, Monoject and Argyle the Cardinal Health LOGO, are trademarks or registered trademarks of Cardinal Health. All other marks are the property of their respective owners. CA Lit. No. 2PATR18-748943, 19-054 (10/2019)



Cardinal Health Canada
1000 Tesma Way
Vaughan, ON L4K 5R8

shop.cardinalhealth.ca